DENTSPLY

ProMix™

AMALGAMATOR

MODEL 400

MODEL 402 CE

***Operations Manual
Directions for Use***


Figure 1


Figure 2

TABLE OF CONTENTS

| | |
|---|--------------------|
| Figures | Inside Front Cover |
| Technical Details | 2 |
| Glossary of Symbols | 3 |
| Definition of Terms. | 3 |
| Operations Guide | 4 |
| Indications for Use | 5 |
| Contraindications and Warnings | 5 |
| Precautions | 5 |
| Installation. | 6 |
| Operation | 6 |
| Setting Memory Locations | 8 |
| Setting Mixing Time | 8 |
| Installing Capsules | 9 |
| Starting the Mixing Cycle | 9 |
| Stopping/Interrupting the Mixing Cycle. | 9 |
| Removing Capsules | 10 |
| Cleaning and Disinfection | 10 |
| Warranty | 10 |

TECHNICAL DETAILS:

| | | |
|---|--|--|
| AC INPUT RATING: | Model 400 Model 402 | 120V, 50/60Hz, 2.0A 230V, 50/60Hz, 1.2A |
| FUSE RATING: | Model 400 Model 402 | F2.0A/250V T1.6A/250V |
| MOTOR RATING: | 115VAC, 60Hz | |
| PRESET MOTOR SPEEDS: | Low (Turtle) High (Rabbit) | 3700 ± 200 rpm 4200 ± 200 rpm |
| IEC EQUIPMENT CLASS: | I | |
| INSTALLATION CATEGORY (OVER VOLTAGE CATEGORY): | II | |
| POLLUTION DEGREE: | 2 | |
| OPERATION MODE: | Intermittent Operation (30 sec ON/120 sec OFF) | |

OPERATING ENVIRONMENT

NOTE: Equipment not suitable for use in the presence of flammable anesthetic mixture with air or with nitrous oxide.

| | |
|------------------------|------------------------------------|
| Operating Temperature: | 10 to 40 degrees C (50°F to 104°F) |
| Relative Humidity: | 30 to 75% (non-condensing) |
| Atmospheric Pressure: | 525 to 795 mmHg |

TRANSPORT: STORAGE CONDITIONS

| | |
|-----------------------|--------------------------------------|
| Temperature: | -10 to +60 degrees C (22°F to 140°F) |
| Relative humidity: | 10% to 100% (non-condensing) |
| Atmospheric Pressure: | 375 to 795 mmHg |

DIMENSIONS:

| | |
|--------------------|-----------------------|
| Height: | 5 inches (12.7 cm) |
| Width: | 9 inches (22.86 cm) |
| Length: | 9.5 inches (24.13 cm) |
| Weight: | 7 pounds (3.18 Kg) |
| Power Cord Length: | 117.6 inches (300cm) |

GLOSSARY OF SYMBOLS

Attention, Consult Accompanying Documents

Alternating Current

Protective Earth

OFF (Power: Disconnection from the Main Power Supply)

ON (Power: Connection to the Main Power Supply)

START /STOP: Initiates and deactivates the amalgamation cycle

HIGH: rabbit

LOW: turtle

DEFINITION OF TERMS:

Initialization: The function performed by the software of the unit after power up or master reset

Stand-By: The condition when the unit is not performing the amalgamation cycle but has completed the “initialization” period

Run: The condition when the unit is performing the amalgamation cycle

OPERATIONS GUIDE

ProMix™ Amalgamator

INTRODUCTION:

The microprocessor controlled ProMix™ Amalgamator is one of the most advanced amalgam and encapsulated material mixer systems presently available in the dental marketplace.

The ProMix™ Amalgamation Unit provides:

- a) **A Constant Mixing Speed** which is automatically maintained during the complete run/mix cycle.
- b) **Two Preset Speeds** with an LED indication of the selected Low speed (3700 rpm) or High speed (4200 rpm).
- c) **Two Memory Locations (M1 and M2)** which can be programmed for time and mixing speed for your two most commonly used encapsulated materials (to the manufacturer's specifications).
- d) **A Timer** to set the run cycle to a maximum of 30 seconds in 1 second increments for either selected speed.
- e) **A Digital LCD Display** for convenient viewing of timing functions.
- f) **A Unique Capsule Holder** (arm assembly) which will accommodate a large variety of today's commonly used capsules (including those with nozzles).
- g) **A Safety Switch** which automatically turns OFF the run cycle should the ProMix™ Amalgamator cover accidentally be opened.
- h) **Quiet** and smooth running operation.
- i) **A Thermally Protected Motor** to prevent any possible damage caused by overheating.
- j) **Manual Cycle Interrupt Capability** to reduce material waste should the timer inadvertently be set longer than required.

INDICATIONS:

For the trituration of dental alloys, capsule delivered temporary filling materials and cements.

CONTRAINDICATIONS:

None known

WARNINGS:

1. We recommend that the ProMix™ Amalgamator not be used in cases where a patient or user of the ProMix™ Amalgamator has been fitted with an implanted heart pacemaker and has been cautioned against the use of small electrical appliances (such as electric shavers, hair dryers, etc.)
2. The ProMix™ Amalgamator is not suitable for use in the presence of a flammable anesthetic mixture with air or nitrous oxide.

PRECAUTIONS:

1. This product is intended to be used as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
2. Operate the ProMix™ Amalgamator only on a hard, flat surface away from countertop edges. Unsafe to operate in any other manner than specified by DENTSPLY.
3. As with all dental procedures, use standard personal protective equipment (i.e., wear face mask, eyewear or face shield, gloves and protective gown).
4. Do not place the unit on or next to a radiator or other heat source. Excessive heat may damage the unit's electronics.
5. Do not operate the unit in a tightly confined space or corner. Place the unit where air is free to circulate on all sides and beneath it.
6. The translucent acrylic arm cover/door of the ProMix™ Amalgamator has a safety interlock switch. Do not attempt to override the switch and/or operate the ProMix™ with this door open.
7. Should the ProMix™ Amalgamator get damaged in any way, do not attempt to operate the unit unless thoroughly inspected and serviced by qualified DENTSPLY Service personnel.

8. Always use the unit in conjunction with a properly grounded power source. Use only the supplied power cord. Connect only to a power source for which this unit is rated.
9. DO NOT IMMERSE UNIT, POWER CORD OR PLUG IN WATER OR DISINFECTANT SOLUTIONS.
See cleaning and disinfection instructions on page (10).
10. DECLARATION OF CONFORMITY - EU
402 Models conform with the Laboratory Equipment European Directive, 73/23/EEC including amendments.

ADVERSE REACTIONS:

Electrical appliances can interfere with some pacemakers. If the user or patient has been fitted with an implanted heart pacemaker, see Warnings.

STEP-BY-STEP INSTRUCTIONS:

INSTALLATION:

NOTE: Upon receiving the ProMix™ Amalgamator system check the packaging and parts for any damage that may have occurred in transit. If damage is apparent, please contact your dealer immediately.

UNIT POWER-UP (INITIALIZATION):

Plug the supplied detachable power cord into the back of the unit and then into the main power supply.

NOTE: Always use the unit in conjunction with a properly grounded power source.

Activate the Power Switch on the back of the unit.

Upon AC power up, the unit will always be initialized as follows:

- the LCD display will read “00”,
- the motor will be OFF,
- the “LOW” motor speed will be selected and the “LOW” LED indicator will be ON.

If these three conditions are not apparent, do not use the unit. Contact qualified DENTSPLY Service personnel.

NOTE: AC power interrupts and degradation causing improper operating voltages to the microprocessor will always result in a microprocessor master reset and re-initialization of the unit.

OPERATION:

The operation of the ProMix™ Amalgamator can be controlled by the LCD display and the seven (7) push buttons/pads located on the front panel (see Figure 1, inside Front Cover).

"LCD" display

Indicates the set (desired) amalgamation cycle/run time and counts down to "00" during operation.

"Up" arrow

Increases the programmed run time by seconds (to a maximum of 30 seconds).

"Down" arrow

Decreases the programmed run time by seconds.

"High" (rabbit)

Selects the high amalgamation speed (4200 rpm) and the associated LED indicator.

"Low" (turtle)

Selects the low amalgamation speed (3700 rpm) and the associated LED indicator.

"Start/Stop" (I/O)

Starts and stops the amalgamation cycle.

NOTE: The two (2) "memory" locations are set at the factory to commonly used settings. To change the settings to suit individual needs see "Setting Memory Locations" (page 8).

"M1" (Memory #1) Stores and recalls #1 preprogrammed amalgamation speed and time, even in the "power OFF" condition.

NOTE: The Memory #1 location is factory set to 15 seconds at Low speed.

"M2" (Memory #2) Stores and recalls #2 preprogrammed amalgamation speed and time, even in the "power OFF" condition.

NOTE: The Memory #2 location is factory set to 20 seconds at High speed.

SETTING MEMORY LOCATIONS:

The ProMix™ Amalgamator memory locations (M1 & M2) can be programmed to suit the two most commonly used amalgamation cycle times and speed settings of the individual user (**corresponding to the manufacturer's specification of the material/amalgam being mixed**).

- a) Select the desired speed by pressing the appropriate key pad (verify the associated LED illuminates).
- b) Select the desired run time by depressing and holding the up or down arrow button/pad (to scroll through time) until the desired time is displayed on the LCD.
- c) Store the setting by depressing and holding the desired memory location (M1 or M2) button/pad for longer than 3 seconds (the LCD will display "00" for 0.5 second before the desired memory setting will reappear indicating the end of the storage cycle).

SETTING THE MIXING TIME:

NOTE: Refer to the encapsulated material manufacturer's specification to determine the proper mixing time and speed.

If the material/amalgam required mixing time and speed is the same as one of the two (2) stored memory locations (M1 or M2) simply depress that desired Memory key pad and proceed to "Installing Capsules".

If not, proceed as follows:

- a) Select the recommended speed by pressing the appropriate speed button/pad:
 - i) High speed (rabbit)
 - ii) Low speed (turtle)
- b) Select the recommended mixing time (to a maximum 30 seconds) by:
 - i) pressing the appropriate Up or Down arrow button/pad until the recommended mixing time is displayed on the LCD; or holding down the Up or Down arrow button/pad automatically scrolls through the mixing time in one (1) second intervals.

INSTALLING CAPSULES:

The ProMix™ Amalgamator capsule holder is designed for easy use and to automatically retain the capsule when installed properly (refer to Figure 2, inside front cover).

- a) Lift the ProMix™ Amalgamator cover to gain access to capsule holder.
- b) Open the capsule holder by pressing the Arm Button with your thumb while holding the Arm Body Tab with your fingers to keep the assembly stable.
- c) Activate the capsule (if necessary, per material manufacturer's instructions).
- d) Place the capsule into the opened jaws of the Arm assembly.

NOTE: Capsules with nozzles should be placed to the left (such that the nozzle fits in the groove of the capsule holder).

To assure consistent mixes, capsules of the same brand should be oriented in the same position (e.g. top to left) each time.

- e) Release the Arm Button, which automatically retains the capsule.
- f) Close the ProMix™ Amalgamator cover.

STARTING THE MIXING CYCLE/UNIT OPERATION:

- a) To start the mixing cycle, simply press the Start/Stop button.
- b) The LCD display will count down to “00” from the selected run time in one (1) second increments.
- c) The ProMix™ Amalgamator will automatically Stop at the end of the run cycle.
- d) The LCD display will return to the previously selected run time.

STOPPING/INTERRUPTING THE CYCLE:

NOTE: Interrupting/halting the mixing cycle in any of the following manners may make the material being processed unusable (refer to the material/ amalgam manufacturers specification).

To stop the mixing cycle (during processing) before the selected run time is complete:

- a) Press any key pad.
- b) OPEN or lift the ProMix™ Amalgamator cover which activates the units safety device.

The reading on the LCD will return to the previously selected value.
The ProMix™ Amalgamator unit will enter/return to the Stand-By mode.

REMOVING CAPSULES:

- a) Lift the ProMix™ Amalgamator cover to gain access to capsule holder.
- b) Open the capsule holder by pressing the Arm Button with your thumb while holding the Arm Body Tab with your fingers to keep the assembly stable.
- c) Remove the processed capsule.

DISINFECTION:

NOTE: Disconnect the power cord from the ProMix™ Amalgamator unit and power source prior to cleaning.

Prior to disinfection, thoroughly wipe all surfaces and power cord clean with a damp towel (with mild detergent if needed). Properly discard used towel. EPA registered surface disinfectants are acceptable. Water-based disinfectant solutions are preferred. Some disinfectants may cause damage to plastic. The disinfectant manufacturer's directions should be followed for optimum results. Do not allow water, cleaning or disinfectant solutions to get inside the unit.

WARRANTY:

The ProMix™ Amalgamator is designed for use in a dental office only. The warranty does not cover other uses. The ProMix™ Amalgamator is warranted against defective material and workmanship for a period of twelve (12) months from the date of purchase by a dentist. In the event of such a defect, DENTSPLY will repair or replace the ProMix™ Amalgamator or necessary parts therein, at its discretion, and such repair or replacement shall be the sole remedy of this warranty. This warranty extends only to the original purchase and is subject to these conditions.

1. Any servicing of the ProMix™ Amalgamator must be performed by qualified DENTSPLY service personnel.

2. The ProMix™ Amalgamator must not be subjected to abuse or improper installation or application, nor installation or repair service by other than trained DENTSPLY service personnel.

THERE ARE NO WARRANTIES, EXPRESS OR LIMITED, WHICH EXTEND BEYOND THIS DESCRIPTION. DENTSPLY neither assumes, nor authorizes any person to assume for it, any other liability in connection with the sale and use of the ProMix™ Amalgamator.

DAMAGES ARE LIMITED STRICTLY TO REPLACEMENT OF THE ProMix™ Amalgamator. DENTSPLY EXPRESSLY DISCLAIMS LIABILITY FOR INCIDENTAL AND CONSEQUENTIAL DAMAGES RESULTING FROM THE USE OF THE ProMix™ Amalgamator.

Claims covered by this warranty will be honored when this warranty is presented through your DENTSPLY Service Center (see list on page 12) within one (1) month from discovery of defect.

Caution: U.S. federal law restricts this device to sale by or on the order of a dentist.

Worldwide Service Centers

DENTSPLY Professional
Attn.: Product Service
1301 Smile Way
P.O. Box 7807
York, PA 17404-0807, USA
1-800-989-8826

DENTSPLY DeTrey GmbH
De-Trey-Str. 1
78467 Konstanz
Germany
Tel. (49) (07531) 583-0

DENTSPLY Limited
Hamm Moor Lane
Addlestone
Weybridge, Surrey KT15 2SE
Tel. (44) (01932) 853422

DENTSPLY DeTrey
ZA du pas du Lac
17, rue Faraday
F-78180 Montigny le Bretonneux, France
Tel. (33) (01) 30147777

DENTSPLY Italia s.r.l.
Via A. Cavaglien 26
1-00173 Roma, Italy
Tel. (39) (06) 7233626

Swiss Representative
DENTSPLY DeTrey
Baar Office
Oberdorfstrasse 11
6324 Baar
Switzerland
Tel. (41) (041) 7662066

DENTSPLY Australia Pty Ltd
11-21 Gilby Road
Mount Waverley, Victoria 3149
Australia

DENTSPLY Sankin K.K.
14-9, Yushima 3-chome
Bunkyo-ku, Tokyo 113-0034
Japan

Hong Kong Division of
DENTSPLY International Inc.
23/F Gee Ghang Hong Center
65 Wong Chuk Hand Road
Aberdeen, Hong Kong

DENTSPLY Industria e Comercio Ltda.
Rua Alice Herve 86
Bengen
25665-010 Petropolis, R.J., Brazil

Dentsply Mexico, S.A. de C.V.
Calzada Vallejo No. 846
Col. Industrial Vallejo
02300 Mexico, D.F.

Amalgam Mixing Record

| ITEM | AMALGAM | MIXING TIME | MIXING SPEED |
|------|---------|-------------|--------------|
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

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